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
EN KERARING - INTRASTROMAL CORNEAL RING SEGMENTS

REV 02

PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

WARNING! This device must only be implanted by an experienced eye surgeon trained in the implantation technique of this specific device.

DESCRIPTION

Keraring – intrastromal corneal ring segments (device, implant) is a medical device implanted in the human corneal stroma for the reduction of corneal surface irregularities and refractive errors associated with corneal ectatic disorders. The device is manufactured from polymethylmethacrylate (PMMA), and is implanted in the corneal stroma as per specific surgical technique. The device acts upon the corneal tissue by altering its central curvature and shape, thus reducing or eliminating morphological irregularities and existing myopia and astigmatism. The device is composed of one semi circular segment of variable arc lengths (90°, 120°, 160°, 210°), variable apical diameters (5mm, 6mm), variable thickness (150µm, 200µm, 250µm, 300µm, 350µm), having a fixed triangular cross section of 600 microns base. Each ring segment tip has a 0.2mm diameter hole to facilitate manipulation and implantation.

INDICATIONS

The correct indication for implantation of the Keraring requires a thorough evaluation of topographic and pachymetric conditions of the cornea, in addition to complete ophthalmologic examination. In general, implantation of the device may be considered in the following cases:

- Keratoconus in contact lens intolerant patients;

- Progressing keratoconus;

- Irregular astigmatism following penetrating keratoplasty;

- Corneal ectasia following previous photoablative procedures (Lasik, PRK, Lasek, Epi-Lasik);

- Irregular astigmatism following radial keratotomy;

- Pellucid marginal degeneration;

- Corneal surface irregularities post trauma resulting in irregular astigmatism

CONTRA-INDICATIONS

Implantation of Keraring must not be considered in the following conditions:

• The safety and effectiveness of Keraring implantation have been established when used as directed. However, visual rehabilitation may be unsatisfactory if this device is not used in accordance to guidelines herein listed, or if patient fails to receive proper post-operative care.

• Some patients with large dilated pupil diameters are predisposed to low light visual symptoms and should be appropriately advised.

• It is anticipated that spectacle correction or contact lens correction may still be necessary following implantation of this device, especially in moderate to advanced keratoconus patients.

• Under mesopic conditions, patients may experience some loss in contrast sensitivity at low spatial frequencies.

• The safety and effectiveness of alternative refractive procedures following implantation of the Keraring has not been established.

HOW SUPPLIED

Each Keraring contains one single intracorneal ring segment, presented sterilized by ethylene oxide in a plastic case sealed inside a sterilization pouch. Device identification labels are provided and indicate model number, device dimensions, implant serial number, and sterilization expiration date. One Keraring patient identification card and is also supplied.

PRECAUTIONS FOR USE AND STORAGE

• Do not use after sterilization expiration date indicated on package.

• Single use only. Do not resterilize.

• This device must only be implanted by an experienced corneal surgeon trained in intrastromal corneal ring implantation techniques and management of related complications.

• Manipulating this device requires the use of the appropriate techniques and instrumentation. Any device

damaged during handling should not be implanted.

• Keraring implantation requires the use of special surgical instrumentation available from Mediphacos. Do not implant or try to implant this device with any other instrument.

• The correct selection of device dimensions and position for a specific patient depends on preoperative manifest refraction and the presence and extent of other corneal conditions, such as keratoconus. Refer to clinical guidelines and updated nomograms available from Mediphacos. Please contact keraring@mediphacos.com

• Do not use if sterile pouch is found opened and/or damaged.

• The device dimensions and sterilization expiration date should be verified before opening the sterile pouch.

• Store at room temperature. Avoid high temperatures above 45° Celsius.

• If for any reason the device is removed from sterile condition and not implanted, the device must not be used and appropriately discarded on hospital waste.

• Implant identification labels supplied must be affixed to patient and/or hospital and/or surgeon files in order to allow tracing each individual device to its end user (patient) if necessary, as required by national and international regulations.

• Implantation of this device must be executed in sterile operating room will all precautions associated with intraocular surgical procedures.

PATIENT PREPARATION AND INSTRUCTIONS

• If patients wears contact lenses, he/she should be advised to discontinue its use at least 2 weeks before their final preoperative examination in order to allow accurate refraction and corneal topography measurements;

• Patients should be instructed not wear any make up, especially around the eyes, on the operative day;

• Patients should be instructed not to rub their eyes for at least 6 months after Keraring implantation, as to avoid implant migration or extrusion and ensure proper incision healing;

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WARRANTY AND LIMITATION OF LIABILITY

Mediphacos Ltda. warrants that the product when delivered is free from defect in materials and workmanship and is in conformity to the manufacturer's current version of product's specification. This warranty is valid until the product expiration date indicated on the external labeling. At its sole option, Mediphacos will replace or provide a refund for any product found to be defective, as long as the defective product is returned to Mediphacos. Mediphacos shall not be liable for any incidental or consequential loss, damage or expense resulting directly or indirectly from the use or the inability to use this product in any given patient. Mediphacos Ltda. neither assumes, nor authorizes any person to assume on its name, any other liability or responsibility with respect to this product other than set forth herein.

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